

## SPECIFICATIONS

Product Description: **FIT TEST KIT**  
Part Number: **2040**  
Style: **SWEET (SACCHARIN)**

### GENERAL DESCRIPTION:

This kit offers qualitative OSHA compliant fit testing for disposable and reusable dust and mist respirators. Saccharin has a sweet taste to indicate failure on the respirator fit test. New disposable glass ampules make pouring solution into nebulizers easier and reduce the chance of contamination.

### KIT:

- Test hood
- Sensitivity nebulizer
- Test nebulizer w/ atomizer
- Wire cleaner
- Box of Sensitivity Solution, 6 ampules
- Box of Test Solution, 6 ampules

### HOOD:

- Poly-coated
- Clear viewing window
- Surged seams
- Reusable

### SWEET (SACCHARIN):

- Glass ampules
- With pre-drilled plastic cap for easy pouring
- Plastic protection from broken glass
- One-time use
- Averages 3-5 fit tests per tube
- Contains 2.5cc of solution

### NEBULIZER:

- Buna-N (Nitrile rubber) air bulb for easy operation
- Converts 5cc (ml) of liquid into aerosol 0.3 to 5.0 microns particle size
- Complete with extra jet and replaceable plastic caps
- Break-resistant polycarbonate plastic (Lexan<sup>®</sup>) that easily disassembles for washing
- Dimensions: 2" x 2" x 7" (5.0 cm x 5.0 cm x 17.7 cm)

### AVAILABLE:

**2040-11K** Sweet (Saccharin) Sensitivity Solution Ampules (6/box)  
**2040-12K** Sweet (Saccharin) Test solution Ampules (6/box)  
**2040-08N** Nebulizers (each)

